

### DNA Isolation Lab—Strawberry

#### 1. Procedures for isolation

- Obtain a Strawberry.
- Place the strawberry in a ziplock bag and add 10 mL of 0.9% NaCl.
- Pulverize the strawberry, in the bag, for about 2 minutes.
- Add 1 mL of detergent such as *Lemon Fresh Joy, Dawn or Woolite*.
- Continue pulverizing for 2 minutes.
- Filter the solutions through several layers of cheesecloth into a clean 50 mL glass beaker (It may be necessary at this point to add an additional 5-10 mL of 0.9% NaCl in order to have enough liquid to filter).
- Tilt the beaker and gently add an equal volume of ice-cold ethanol down the inside of the beaker.
- A stringy, white precipitate containing DNA will form along the interface between the two solutions.
- Place a wooden skewer upright in the center of the beaker and gently spool the DNA by gently turning the skewer continuously in the same directions.
- The collected DNA sample may be stored in a small container with a small amount of the 0.9% NaCl solution.

#### Discussion Questions:

1. What organelle originally held the sticky white precipitate?
2. If you had “molecular vision” and could see the molecules that make up the stingy white precipitate,
  - a. what elements would you see in it? \_\_\_\_\_
  - b. what subunit would you see in it? \_\_\_\_\_
  - c. what sugar would you see? \_\_\_\_\_
  - d. what possible bases could you see? \_\_\_\_\_
3. How is the structure of DNA related to its function?
4. Through what process do cells like those found in plants make more DNA?

**BONUS (10 pts – Use the AP Biology Textbooks!)) Label the replication fork shown below:**

